

Deluxe Viral Disinfection Kit



Inactivates Viral Pathogens

This kit is designed to allow anyone to quickly make disinfectant solutions which are effective against contagious viral and bacterial pathogens including Human Corona Virus, Influenza, Adenovirus and many others.

64 Millennium Q (EPA Reg no. 1839-95-5449) has demonstrated effectiveness against viruses similar to 2019 novel coronavirus (SARS-CoV-2) on hard, non-porous surfaces. Therefore, this product can be used against SARS-CoV-2 when used in accordance with the directions for use against Norovirus (Norwalk Virus) on hard, non-porous surfaces.

Check List of High Touch Critical Surfaces:

- Bedrails
- Tables
- Light switches
- Door knobs
- Hand rails
- Window sills
- Desk tops
- Chairs
- Telephones
- Door frames



Makes 4 full quarts of Hospital Grade Disinfectant.

The Kit includes:

- 1 Viral disinfectant concentrate
- 1 Mixing Bottle
- 4 Empty quart bottles
- 4 Microfiber Towels
- Detailed Instructions



Fill mixing bottle with water up to within 2.5 inches from top. Empty the 2 oz bottle of Millennium Q into the mixing bottle, place cap on bottle, gently mix.



Fill quart bottles with ready to use solution you just made. Attach trigger sprayer.



Spray all hard surfaces, allow to remain wet for ten minutes.