

Deluxe Viral Disinfection Kit

Inactivates Viral Pathogens Kit Part number 903020



This kit is designed to allow anyone to quickly make disinfectant solutions which are effective against contagious viral and bacterial pathogens including Human Corona Virus, Influenza, Adenovirus and Norovirus.

64 Millennium Q (EPA Reg no. 1839-95-5449) has demonstrated effectiveness against viruses similar to 2019 novel coronavirus (SARS-CoV-2) on hard, non-porous surfaces. Therefore, this product can be used against SARS-CoV-2 when used in accordance with the directions for use against Norovirus (Norwalk Virus) on hard, non-porous surfaces.



Apply RTU label to mixing bottle. Fill mixing bottle with water to within 2.5 inches from the top as shown.



Empty the 2 oz bottle of Millennium Q into the mixing bottle, cap and gently agitate bottle.



Fill quart bottles with the ready to use solution you just made.



Attach trigger sprayer or flip top cap.



Apply directly or pre-moisten cloth and wipe surface. Surface must stay moist for 10 minutes. Allow to air dry.



Wipe bedrails, tables, light switches, door knobs, hand rails, window sills, desk tops, chairs, phones, and any other high touch surfaces.

Multi-Clean®

600 Cardigan Road
Shoreview, Minnesota
Phone 651-481-1900
Web www.multi-clean.com

See 64 Millennium Q SDS
for precautionary statements.

004204 08042020